

Uterine Contour™


- Uterine Contour™ automatically detects endometrium and provides a clear coronal view of the uterus
- It provides 2 classifications (ESHRE/ESGE, ASRM) to analyze the uterine shape
- The classification of uterine malformation is assigned to the report



How to enter the mode

Method 1

1. Acquire the sagittal uterine view with 3D endovaginal probe
2. Press [Uterine Contour] button
3. Contour ROI is automatically drawn and 3D coronal view is displayed



Method 2


1. Acquire the sagittal uterine view with 3D endovaginal probe
2. Acquire 3D volume data
3. Touch the [Oblique] tab
4. Touch the [Auto Contour] button
5. Contour ROI is automatically drawn and 3D coronal view is displayed

* Uterine Contour is available in Adnexa, General and General1 presets of GYN app.



Further manipulation

1. Modify the contour line with [Rotate Line] and [Line Position]
 - [Rotate Line] : Adjust the contour line by rotating the knob key
 - [Line Position] : Activated by pressing the [Change] key
2. [OVIX] shows the volume data with thickness information
3. Press [Classification] and select the uterine shape on touch panel to assign on the report
 - Report page after [Classification]



Ultrasound Report			
Name	ID	Exam. Date	2020-12-14
[Gynecology]			
Abnormalities of the uterus		Value	
Classification (ASRM)		I. Hypoplasia /agenesis - Cervical	



Setup Guide

1. Classification guideline can be set in Measure Setup > Measure > Gynecology
2. There are two types of Uterine Classification
 - ESHRE/ESGE
 - ASRM

