

S-Shearwave Imaging™

Calculates speed of transvers wave inside of the ROI and displays stiffness index by color coded map. It increases the accuracy of clinical diagnose such as diffuse liver disease.

● Preparation

- ▶ Fast overnight before the examination. (At least 4 hours)


● Obtaining B-mode

- ▶ Acquire plane perpendicular to the right lobe capsule in supine or slight left lateral position with right arm extension.
- ▶ Contact a transducer well to intercostal space to get a high SWI SNR image.
- ▶ Scan the right intercostal space avoid shadowing area and adjust Gain not to be too bright or too dark image.

● Stiffness map Positioning


- ▶ The Map should be positioned at least 1.5cm-2cm below the liver capsule and avoid reverberation artifact. *(Depth of ROI should be at least twice the thickness of the muscle and fat layer in obese patient)*
- ▶ Stiffness map should be close to the center of the image.
- ▶ Position the stiffness map less than 6cm-7cm depth, it is limitation to measure an adequate shearwave.
- ▶ Avoid visibly large vessels if possible. *(Although the algorithm automatically excludes vessels.)*

Stiffness map Positioning example



[Bad example]

- ① Should avoid reverberation artifact.
- ② Should avoid shadow artifact.



③ Dual Live mode (RMI & Stiffness map)

[Good example]

- ① ROI should be positioned at least 1.5-2cm below the liver capsule.
- ② Stiffness map should be placed less than 6-7cm depth.
- ③ Display [RMI] and [Stiffness] map on 2D image in real time.


● Use RMI(Reliability Measurement Index) Map

- Recommend to display RMI map and Stiffness map together by pressing Dual Live for high reproducibility.
- High RMI values are strongly correlated with reproducible measurements.
Maximum RMI is 1, and it is considered to be reliable RMI 0.4 at least.

● Measurement & Quality criteria

- When if stiffness map is unstable, keep acquiring images for 5 - 6sec with gentle breath.
- At least of 1cm sized measure ROI is recommended.
- The optimal measure ROI location for acquisition depth is 4 - 4.5cm.
- Position a measure ROI in an area that meets both conditions : High RMI area and Good uniformity on stiffness map.
- At least 5 times of measurement are recommended with a quality assessment.
- IQR/Med(kPa) should $\leq 30\%$ or IQR/Med(m/s) $\leq 15\%$.
- Each result requires to have RMI ≥ 0.4 .


Dual Live Mode



[Bad example]

Low reproducibility (Measure on less uniformity area)

- RMI = 0.9, Mean = 9.6kPa



[Good example]

High reproducibility (Measure on good uniformity area)

- RMI = 0.9, Mean = 5.1kPa