

UterineContour™

A feature to extract the centerline and thickness of endometrium semi-automatically

Tips & Tricks

Clinical Values

I Clinical Benefits & Cases

- UterineContour™ can show a coronal view in easy way, which is useful to figure out the female genitalia anomaly.
- Save your manipulation time to find a ideal coronal view to assess an endometrium shape in 3D.
- In addition, uterine malformation classification are reported according to the *ESHRE/ESGE or **ASRM guideline selection.



* ASRM : American Society for Reproductive Medicine
** ESHRE/ESGE: European Society of Human Reproduction and Embryology/ European Society for Gynaecological Endoscopy

Configuration

I Checklist

- System: 100% (100% / 100%)
- System: 100%
- Product: 100% (100% / 100%)

Tips & Ticks

01. Pre-Condition for UterineContour™

- Uterine contour™ works regardless an endometrium thickness, which is changing depends on a reproduction cycle. There are in Early proliferative, secretory, premenstrual and additionally post-menopause phase.
- Recommend to demonstrate endometrium's echogenic centerline in Sagittal view for high detection rate and accuracy.



02. How to do the UterineContour™

- Take a Volumetric Uterus Coronal view through Uterine Contour by Semi-automatic operation.



① Take a sagittal view with an endometrium line.

- ② Tap the [Uterine Contour] button on the Touchscreen.



③ Both A plane and Oblique(OVIX) view are shown.

03. UterineContour™ Adjustment



- Line Position**
 Update a contour line by moving a trackball.
- Chroma**
 Switch a 2D or 3D Plane's Chroma map. Sepia map (Map 2 - 7) is usually preferred for contrasty image.

- Rotate Line** : Rotate a contour line with a knob button.




- Line position and Rotation Line** adjustment would be helpful to navigate and define an uterine congenital anomaly.

04. Measurement

- Measurement assessment is supported with caliper to classify a type of uterine structure : Bicornuate, Septate, Arcuate uterus

- 1) Measure a distance between the internal line and the uterine fundus (blue image)
- 2) Measure the outer surface. Take distance between the external line and present with normal the apex of fundus external contour




05. Assign a Classification

- ESHRE/ESGE Classification
- ASRM(1998) Classification
- ASRM MAC 2021 Classification
- Select your default classification
→ Utility > MeasureSetup > Gynecology > Uterus
- Report Assign by tapping the [Classification]

Type of uterus is assigned to the report description by tapping a specific one on the Touchscreen.



- Report page : The assigned Classification will be shown on the Report



Disclaimer

- Before using the material, proceed the review of the local advertising regulation through your RA and legal team.
- The products, features, options, and transducers may not be commercially available in some countries.
- Sales and shipments are effective only after the approval by the regulatory affairs. Please contact your local sales representative for further details.