



## S-Shearwave Imaging™

### RS85 Prestige Quick Guide



# 1. Start S-Shearwave Imaging™




## 1 S-Shearwave Imaging

Tap the [S-Shearwave Imaging] on the touch screen to start.

## 2 Dual Live

It displays both stiffness map and RMI map or 2D images in real time.




## 3 HQ Mode

With this [HQ Mode], it provides high resolution image.

## 4 Elasticity Range

It is to change the color stiffness map range. If you change the [Elasticity Range], color map within ROI and the color of elasticity bar also will be changed accordingly.

## 2. Parameters



[RMI map]

### 1 RMI

It represents the reliability of the measured the shearwave elasticity value.

- The higher RMI index indicates higher reliability.
- It is considered to be reliable for  $RMI \geq 0.4$

### 2 RMI map

It can change the RMI map for user preference.

### 3 RMI Threshold

User can set the [RMI threshold]. then SWI color map will be excluded below the threshold.

### 4 Unit



It is to change the unit of an elasticity value.

- The unit [kPa] is based on Young's modulus.
- The unit [m/s] is for Shearwave speed.

### 5 Site

It is to change the ROI site.

# 3. Measure the Elasticity (1)




## 1 Freeze with Meas./ Ratio Meas.

Starts Elasticity Measurement/Elasticity Ratio Measure when pressing Freeze.

## 2 Next Frame for Meas.

- Auto: Automatically moves to the next frame after confirming the measurement by pressing [Set] button.
- Manual: Moves to the next frame by pressing [Next Frame] after confirming the measurement.





## Notes

The Next Frame direction can be selected between Forwards or Backwards in the Setup page.

*\*Setup (Imaging ► Features ► S-Shearwave Imaging ► Next Frame Direction ► **Forward / Backward**)*

# 3. Measure the Elasticity (2)



## 1 Elasticity Measure

Tap the [Elasticity Measure] button.

- You can adjust ROI size using change key.
- Press the [Set] button to save the value.
- When moving the trackball, it copies additional ROI for next measurement.

## 2 Elasticity Ratio Measure

Tap the [Elasticity Ratio Measure] button.

- Place the first quantification ROI on the **area of interest** within the elasticity ROI and press the [Set] button.
- Place the second quantification ROI on the **reference** area within the elasticity ROI and press the [Set] button.



## 3 Meas. Type

Select the measurement of ROI type between **Circle** or **Trace** type.

### Notes

If Next Frame for Meas. is used, the location of the measurement ROI to the next frame can be set in the Setup page.

\*Setup (Imaging ▶ Features ▶ S-Shearwave Imaging ▶ Next Measurement Cursor Display Position ▶ **Center of ROI / Keep Last Position**)




The stiffness results within quantification ROI are displayed.

#### 1 Result box

- Mean : Mean Elasticity value in ROI.
- SD : Standard deviation in ROI.
- Min : Minimum Elasticity value in ROI.
- Max : Maximum Elasticity value in ROI.
- Diameter : Diameter of ROI box.
- RMI : The reliability of the measured Shearwave Elasticity value.

#### 2 Stiffness Result

Available to view the measured results. It can scroll the result using the [Scroll] button. If you want to delete the result individually by pressing the [Delete] button or selecting  icon.



#### 3 Measure Summary Display

- Available to select if the measure summary (Median, Mean, IQR, IQR/Med value) display on the result or not.
- You can configure the results that you wish to display and their order in 'Utility > Setup > Imaging > Features > S-Shearwave imaging™'.

#### 4 Store with Result

Saves the image including the measurement results area.

## 5.Profile (1)



## 1 Profile


Tap the [Profile] button on the touch screen.

## 2 Results


- S-Shearwave Imaging™ Profile view could be used criterion to assess the quality of the data.
- **Median** : Central value
- **IQR/Med.**
  - IQR(Interquartile Range) : 75percentile - 25percentile
  - IQR/Med < 30% : Acceptable value
  - IQR/Med > 30% : Unreliable value

## 3 Result Graph

Average of shear wave measurements



Standard deviation of measurements



## 1 Auto Profile

Exclude all unreliable data to set the IQR/Med value to 30% or lower.

## 2 Semi-Auto Profile

Exclude unreliable data one by one, start from the most unreliable data, and set the IQR/Med value close to 30% or lower.

## 3 Save

Save and send the result to the report page.

## 4 Report

Select the [Report] button on the touch screen and the result value is saved on the report.

4

Ultrasound Report

Name: SAMSUNG    ID:    Exam. Date: 2020-06-22

[ Abdomen ]

S-Shearwave Imaging

Elasticity [Site A]


ROI	Mean		SD		Min-Max		Diameter cm	RMI
	kPa	m/s	kPa	m/s	kPa	m/s		
1	5.5	1.35	1.2	0.14	3.9 - 9.7	1.14 - 1.80	1.0	0.9
2	5.3	1.33	1.2	0.14	3.5 - 9.3	1.07 - 1.76	1.0	0.8
3	5.2	1.31	0.8	0.10	4.1 - 7.7	1.17 - 1.60	1.0	1.0
4	5.3	1.33	0.8	0.10	3.5 - 7.4	1.07 - 1.57	1.0	0.9
5	5.5	1.35	1.0	0.12	4.1 - 9.7	1.17 - 1.80	1.0	0.9

Summary: Mean: 5.4kPa / 1.33m/s    Median: 5.3kPa / 1.33m/s  
IQR: 0.3kPa / 0.03m/s    IQR/Med.: 4.7% [kPa] / 2.3% [m/s]

## Notes

The result can be trusted only when the value from dividing Interquartile Range by Median (IQR/Med.) is less than 30%.





## Notes

### How to obtain accurate S-Shearwave Imaging™.

- Fast for 4 hours before the examination is recommended.
- **To maximize intercostal space, supine position** with raised a right arm above the head is recommended.
- Scan through the intercostal space between ribs since evaluation on right liver lobe is more recommended.
  - ✂ Possible to get more reliable result than through a subcostal scan by preventing motion artifacts and providing an easy access.
- Ask patients to breathe steadily and stably before holding their breath to avoid deep breathing.
- ROI should be at least 1.5cm below the liver capsule. (over 2cm is recommended)

- The features, options may not be commercially available in some countries.
- Sales and shipments are effective only after the approval by the regulatory affairs. Please contact your local sales representative for further details.
- This Quick guide does not include all of the details of instruction, for more detail, please refer to RS85 User Manual.
- Do not distribute this document to customers unless relevant regulatory and legal affairs officers approve such distribution.
- This User Quick Guide is based on RS85 V2.08.01
- Disclaimer: Some Images in this content were obtained from other system.

**SAMSUNG MEDISON CO., LTD.**

© 2024 Samsung Medison All Rights Reserved.

Samsung Medison reserves the right to modify the design, packaging, specifications, and features shown herein, without prior notice or obligation.