


5D NT™

HERA Z20 Quick Guide



1. Activate 5D NT (Before 3D Acquisition)




1	3D Button	Press [3D] button on the control panel.
2	5D Ready	Tap [5D Ready] on the touch panel.
3	5D NT	Tap [5D NT] to activate 5D NT™.
4	Preset	Select preferred preset between default and users. Each 5D feature could have own preset parameters.
5	Parameters	[Angle] and [Quality] are customizable.
6	Volume acquisition	Find fetal NT on sagittal plane as clear as possible.


★ Tips

For more accurate NT evaluation, follow NT measurement guide-line created by FMF (Fetal Medicine Foundation).


2. Activate 5D NT (After 3D Acquisition)



① 5D	Tab to activate [5D] features.
② 5D NT	Tab [5D NT] button to activate 5D NT™.



3. Set cursor at Diencephalon



1 Seed Point


Set [A] cursor on fetal diencephalon. Then system automatically finds true mid-sagittal plane.

2 Auto MSP

In case of incorrect mid-sagittal plane, turn off [Auto MSP] so that you will set the cursor on Diencephalon once again.

3 Auto Zoom

When it is turned on, it automatically magnifies the image.





★ Tips

Slices of a 3D image will be displayed as images rendered in 3D. You may select [Slab 3D Type] among Surface, X-ray, Max and Min.

5D NT™

4. Measure NT or IT



<p>① Angle</p>	<p>Angle difference from initial acquisition.</p> <ul style="list-style-type: none"> ▪ Blue Solid Line : MSP angle ▪ Red Dotted Line : Initially acquired angle
<p>② Caliper Set</p>	<p>Place caliper on fetal NT or IT depending on selection of item.</p>
<p>③ Result</p>	<p>Measurement result at the maximum distance of NT or IT.</p>
<p>④ Item</p>	<p>Select type of item between NT and IT</p>
<p>⑤ Caliper Placement</p>	<p>Select type of caliper placement among "On to On", "On to Max Brightness", and "In to In".</p>
<p>⑥ Assign to Report</p>	<p>Tap to apply calculated result to the report.</p>

- The features, options may not be commercially available in some countries.
- Sales and shipments are effective only after the approval by the regulatory affairs. Please contact your local sales representative for further details.
- This Quick guide does not include all of the details of instruction, for more detail, please refer to HERA Z20 User Manual.
- Do not distribute this document to customers unless relevant regulatory and legal affairs officers approve such distribution.
- This User Quick Guide is based on HERA Z20 V1.00.
- Disclaimer: Some Images in this content were obtained from other system.

SAMSUNG MEDISON CO., LTD.

© 2024, Samsung Medison All Rights Reserved.

Samsung Medison reserves the right to modify the design, packaging, specifications, and features shown herein, without prior notice or obligation.