


RealisticVue™

HERA Z20 Quick Guide



1. Probe and Preset



※ RealisticVue™ can be operated under the following conditions :

Probe	Application
CV1-8A EV2-12 CMV1-10	All presets

2. Volume acquisition



① Volume acquisition

Acquired volume data with the structure you wish to apply RealisticVue.

★ Tips


For fetal face 3D evaluation, obtain fetus profile as much as clear without interruption of hand or cord.

RealisticVue™

3. Apply RealisticVue™




- 1 **RealisticVue** Tap [RealisticVue] to apply it.
- 2 **Light Direction** Tap [Light Direction] within short cut keys when you adjust the direction of light freely by trackball.
- 3 **Render Setup** Tap [Render Setup] on the touch screen to adjust parameters of RealisticVue.



- 4 **Color Set** Select color of RealisticVue effect among 8 designated options.
- 5 **Color parameters** For fine control of color effect, adjust Hue, Saturation and Lightness.
 - **Hue:** Shade of a color
 - **Saturation:** Purity or intensity of color
 - **Lightness:** Brightness of color

3. Apply RealisticVue™



1 Pre-Defined Light Direction


Tap preferred light direction among 9 designated options.

2 Clear SFVI

Smart Filter Volume Imaging. It reduces noise.

RealisticVue™

4. Color Set of RealisticVue™



- The features mentioned in this document may not be commercially available in all countries. Due to regulatory reasons, their future availability cannot be guaranteed.
- Do not distribute this internal document to customers unless relevant regulatory and legal affairs officers approve such distribution.
- This product is a medical device, please read the user manual carefully before use.
- This document is provided to help you understand the feature.
- This User Quick Guide is based on HERA Z20 V1.00.
- Disclaimer: Some Images in this content were obtained from other system.

SAMSUNG MEDISON CO., LTD.

© 2024 Samsung Medison All Rights Reserved.

Samsung Medison reserves the right to modify the design, packaging, specifications, and features shown herein, without prior notice or obligation.