

## WHITE PAPER

# MV-Flow™ and its quantitative Vascular Index (VI) for differential diagnosis of benign or malignant adnexal mass

**Prof. Hong Xiang**  
**Linlin Ruan**

Obstetrics and Gynecology Ultrasound Department,  
First Affiliated Hospital of Xinjiang Medical  
University, Xinjiang Key Laboratory of Ultrasound  
Medicine, Urumqi, China

## Introduction

### • Background

Ovarian cancer ranks seventh among cancers affecting women globally and is the eighth leading cause of cancer-related deaths<sup>1</sup>. Ultrasound(US) is the preferred initial imaging modality for screening ovarian-adnexal tumors<sup>2,3</sup>. Color Doppler can visualize tumor blood flow, and Pulsed wave(PW) Doppler with a resistance index(RI)  $\leq 0.4$  indicates a higher likelihood of malignancy<sup>4</sup>. However, Doppler technology has limitations regarding blood flow velocity and angles. Contrast-enhanced ultrasound (CEUS) offers greater sensitivity and accuracy in distinguishing between benign and malignant ovarian tumors compared to conventional grayscale and Doppler US, albeit with lower specificity<sup>5</sup>. Furthermore, CEUS is costlier, requires extended observation periods, and carries potential risks of adverse reactions. MV-Flow™ is the technique that detects microvascular flow undetectable by conventional Doppler imaging. Unlike color or power Doppler, where low-velocity blood flow can be mistaken for background tissue or noise, MV-Flow™ offers suitable sensitivity and resolution by suppressing tissue noise. Meanwhile, MV-Flow™ offers a semi-quantitative measuring tool, the Vascularity Index(VI<sup>MV</sup>). VI<sup>MV</sup> is a ratio of the MV-Flow™ pixels(or area) to the total pixels(or area) in the region of interest(ROI). It represents the abundance of blood vessels within the ROI. In the absence of a standardized criterion for MV-Flow™ interpretation, VI<sup>MV</sup> can serve as a useful judgment index for evaluating MV-Flow™ findings. Ultimately, MV-Flow™ excels in detecting microvascular flow within tissues and organs, surpassing the capabilities of conventional Doppler technology.

### • Purpose

The purpose of this study was to observe the value of MV-Flow™ for differential diagnosis of ovarian-adnexal masses.

## Methods

### • Study design and population

The study was performed under ethics approval from Institutional Clinical Research Ethics Committee of the First Affiliated Hospital of Xinjiang Medical University. Written informed consent was obtained from all participants.


From June 2021 to July 2022, a total of 107 patients with ovarian-adnexal masses who were treated at this hospital were included. The inclusion criteria are: ① Ovarian-Adnexal Reporting and Data System (O-RADS) diagnosed as category 2 to 5<sup>6</sup>; ② no treatment before the transvaginal US examination; ③ within 7 days before the operation; ④ definitive diagnosis could be confirmed by surgical pathology or by follow-up.

### • Ultrasound parameters acquisition

All patients underwent transvaginal US examination and MV-Flow™ imaging by two US doctors who had more than 8 years practicing experience in obstetrician and gynecologic sonography, using Samsung HERA W10 diagnostic ultrasound system with a transvaginal probe EV3-10B.


The patient lies supine on the examination table. The uterus and adnexa were scanned from multiple planes and angles. When lesions were present in both adnexal regions, if their US appearances were similar, the one with the larger maximum diameter or more complex sonographic features was selected for analysis; if different, each was included separately.

First, routine US examination of the mass was performed and the mass was categorized according to O-RADS. Then MV-Flow™ images were retained from areas with abundant blood flow (Figure 1.a). The US doctors outlined the largest area of the tumor using Manual Trace to obtain its  $VI^{MV}$  value (Figure 1.b). If the tumor is too large and exceeds the probe's scanning range, the area with the most abundant blood flow should be displayed, and the lesion's maximum area should be marked as much as possible. For complex cystic tumors, if the echogenicity of solid component within was  $< 50\%$ , then measure the  $VI^{MV}$  of the solid component (Figure 2.a). If the solid component has an echogenicity  $\geq 50\%$ , then measure the VI of the entire tumor (Figure 2.b). Each US doctors measured each lesion at least 3 times, and selected the maximum  $VI^{MV}$  from the respective measurements as the result.




**Figure 1.** An MV-Flow™ image of a case of high-grade serous carcinoma of the left ovary, for illustrating how we obtained the section and  $VI^{MV}$  values.

a. MV-Flow™ images were retained from areas with abundant blood flow.




b. Outlined the largest area of the tumor using Manual Trace to obtain its  $VI^{MV}$  value.



**Figure 2.** The different methods of drawing the ROI for VI<sup>MV</sup> in tumors with varying proportions of solid components.

a. If the echogenicity of solid component within was < 50%, then measure the VI<sup>MV</sup> of the solid component (this case is a right ovarian serous carcinoma).



b. If the solid component has an echogenicity ≥50%, then measure the VI<sup>MV</sup> of the entire tumor (this case is a left ovarian borderline serous cystadenoma).

### • Pathological Histological Examination

The benign or malignant diagnosis of the tumor was confirmed by surgical pathology or through follow-up.

### • Statistical analysis

MedCalc 19.0 statistical analysis software was used to evaluate the consistency of the measurement of VI<sup>MV</sup> between the 2 US doctors by intra-class correlation coefficient(ICC), and ICC>0.75 was considered to be a high level of consistency, and the data could be used for subsequent analysis. The median (upper and lower quartiles) was used to express the data that did not conform to normal distribution, and non-parametric tests were performed. The receiver operating characteristic (ROC) curve was plotted and the area under the curve (AUC) was calculated to evaluate the value of differential diagnosis of benign and malignant ovarian-adnexal tumors by VI<sup>MV</sup> obtained from MV-Flow<sup>TM</sup> imaging, and a difference of P<0.05 was considered statistically significant.

## Results

### • Participant characteristics

A total of 107 patients, aged between 22 and 81 years, with a mean age of 42.2 ± 13.2 years, were enrolled. Additionally, 115 ovarian-adnexal tumors were included, comprising 33 malignant (including borderline) and 82 benign tumors. One benign lesion disappeared after three follow-up visits, confirmed as enlarged adnexal tissue, while the remaining 114 tumors were verified through surgical pathology (Table 1).

**Table 1.** The pathological results of 114 surgically confirmed

Benign masses (n=81)		NO.	Malignant masses (n=33)		NO.
Endometrial implantation cyst	29		Serous carcinoma	18	
Mature teratoma of ovary	16		Granulose cell tumor	3	
Serous cystadenoma	7		Borderline tumor	3	
Serous Adenofibroma	7		Mucinous carcinoma	2	
Pure ovarian cyst/Ovarian crown cyst	6		Clear cell carcinoma	3	
Functional cysts/Luteal cysts	3		Immature teratoma	2	
Serous cyst	2		Endometrioid carcinoma	1	
Mucinous cystadenoma	2		Dysgerminoma	1	
Fallopian tube abscess	2				
Thecoma	2				
Plasma-mucinous cystadenoma	1				
Laevicellulare	1				
Accessory ovary	1				
Tuberculosis	1				
Broad Ligament Leiomyoma	1				



• The MV-Flow™ imaging results

MV-Flow™ imaging showed that in malignant ovarian-adjacent tumors, blood flow appeared mostly linear and centrally distributed, whereas in benign tumors, blood flow was predominantly punctate or showed limited linear flow with a peripheral distribution.




Figure 3.

a/b. A malignant Granulosa cell tumor with color Doppler & MV-Flow™ showing blood flow centrally distributed.  
c/d. A benign Serous Adenofibroma with color Doppler & MV-Flow™, showing the limited blood flow peripherally distributed.

• **The VI<sup>MV</sup> diagnostic efficacy**


The ICC of the tumor VI<sup>MV</sup> measured by two doctors is 0.987 (P<0.001), indicating high consistency.

The median VI<sup>MV</sup>(%) of malignant ovarian-adenexal tumors was 16.57 (8.85, 24.50), while that of benign tumors was 1.20 (0.58, 3.05). The difference was statistically significant (Z=6.995, P<0.001).

VI<sup>MV</sup>(%) for distinguishing between benign and malignant ovarian-adenexal tumors had an AUC of 0.918 (0.852, 0.961). Using a cutoff value of VI<sup>MV</sup>(%) ≥ 7.15, the accuracy for diagnosing malignant masses was 89.57% (103/115), with a sensitivity of 84.85% (28/33), specificity of 91.46% (75/82), positive predictive value of 80.00% (28/35), negative predictive value of 93.75% (75/80), positive likelihood ratio of 9.94, and negative likelihood ratio of 0.17.

**Table 2.** VI<sup>MV</sup>(%) for Differential Diagnosis of Benign and Malignant Ovarian-adenexal Tumors

VI <sup>MV</sup> (%)	Pathology (NO.)		SUM
	Malignant	Benign	
≥7.15	28	7	35
<7.15	5	75	80
SUM	33	82	115



**Figure 4.** ROC curve for differential diagnosis of benign and malignant ovarian-adenexal tumors using VI<sup>MV</sup>(%)




## Discussion

Using ultrasound to assess the abundance of blood flow signals is one of the important imaging methods for differential diagnosis of benign and malignant ovarian-adnexal tumors, and the O-RADS system also includes content on color scoring<sup>6</sup>.

This study showed that different adnexal tumors exhibit different characteristics on MV-Flow<sup>TM</sup> imaging. Malignant tumors tend to have abundant blood supply, complex and relatively tortuous blood vessels, more branching vessels, and the main blood vessels are often distributed near the center of the mass or within solid components. In contrast, benign tumors generally possess less blood supply, simpler vessel pathways, with blood vessels distributed around the periphery or within the septa of the mass, often appearing as punctate or scanty linear structures.

MV-Flow<sup>TM</sup> imaging-derived vascularity index (VI<sup>MV</sup>) also demonstrated high diagnostic efficacy in distinguishing between benign and malignant adnexal tumors. With a cutoff value of 7.15(%), the AUC for distinguishing benign from malignant tumors was 0.918, with a sensitivity of 84.85% and specificity of 91.46%.

In another similar study by my team, we compared the diagnostic performance of VI<sup>MV</sup> with O-RADS system. Using cutoff values 7.15(%) for VI<sup>MV</sup> and categorizing O-RADS into 5 classes, the AUCs for distinguishing between benign and malignant ovarian-adnexal tumors were 0.923 and 0.929 separately. They were very similar. When these two methods were combined, the AUC was 0.955 (Figure 5)<sup>7</sup>.



**Figure 5.** Receiver operating characteristic curves of O-RADS US and MV-Flow (VI<sup>MV</sup>) alone and in combination for diagnosing benign and malignant ovarian-adnexal masses. O-RADS US, Ovarian Adnexal Reporting and Data System Ultrasound; VI<sup>MV</sup>, Vascular index from MV-Flow<sup>7</sup>.

However, the O-RADS system is complex, and MV-Flow™ and VI<sup>MV</sup> measurements are relatively simple, consistent, and easy for young doctors to master.

In this study, we also found that when used for small tumors (e.g., diameter < 10 cm), MV-Flow™ imaging exhibited fewer artifacts and smaller errors in VI measurements. Thus, MV-Flow™ imaging is more valuable for distinguishing relatively small-volume ovarian-adnexal tumors.

## Conclusion

MV-Flow™ imaging and its quantitative VI<sup>MV</sup> are highly effective in the differential diagnosis of benign and malignant ovarian-adnexal tumors.

## References

1. Jung SI. Ultrasonography of ovarian masses using a pattern recognition approach. *Ultrasonography* 2015;34:173-182.
2. Yoeli-Bik R, Lengyel E, Mills KA, Abramowicz JS. Ovarian masses: the value of acoustic shadowing on ultrasound examination. *J Ultrasound Med* 2023;42:935-945.
3. Vazquez-Manjarrez SE, Rico-Rodriguez OC, Guzman-Martinez N, Espinoza-Cruz V, Lara-Nunez D. Imaging and diagnostic approach of the adnexal mass: what the oncologist should know. *Chin Clin Oncol* 2020;9:69.
4. SEHGAL N. Efficacy of color Doppler ultrasonography in differentiation of ovarian masses[J]. *J Midlife Health*, 2019,10(1):22-28.
5. XUN L, ZHAI L, XU H. Comparison of conventional, Doppler and contrast-enhanced ultrasonography in differential diagnosis of ovarian masses: A systematic review and meta-analysis[J]. *BMJ Open*, 2021,11(12):e052830.
6. ANDREOTTI R F, TIMMERMAN D, STRACHOWSKI L M, et al. O-RADS US risk stratification and management system: A consensus guideline from the ACR ovarian-adnexal reporting and data system committee[J]. *Radiology*, 294(1):168-185.
7. Ruan L, Liu H, Xiang H, Ni Y, Feng Y, Zhou H, et al. Application of O-RADS US combined with MV-Flow to diagnose ovarian-adnexal tumors. *Ultrasonography* 2023 Aug 25 [Epub].<https://doi.org/10.14366/usg.23061>

## Disclaimer

- \* The features mentioned in this document may not be commercially available in all countries. Due to regulatory reasons, their future availability cannot be guaranteed.
- \* Do not distribute this document to customers unless relevant regulatory and legal affairs officers approve such distribution.
- \* Images may have been cropped to better visualize their pathology.
- \* This clinical practice review is a result of a personal study conducted by collaboration between Samsung Medison and Prof. Hong Xiang & Linlin Ruan.
- \* This review is to aid customers in their understanding, but the objectivity is not secured.
- \* 본 자료는 삼성메디슨이 연구자와 협업하여 산출된 개인 연구의 결과물입니다.  
고객의 요청에 따라 이해를 돕기 위해 제공하는 자료일 뿐 객관성은 확보되지 않았습니다.



Scan code or visit  
[samsunghealthcare.com](https://samsunghealthcare.com)  
to learn more

---

**SAMSUNG MEDISON CO., LTD.**

©2024 Samsung Medison All Rights Reserved.

Samsung Medison reserves the right to modify any design, packaging, specifications and features shown herein, without prior notice or obligation